

BIOIMAGING SHARED RESOURCE
POWER ON/OFF PROCEDURES

SYSTEMS

BLSB 910

- Leica SP8 tau-STED*
- DeltaVision Ultra Widefield*

BLSB 912

- Nikon A1R-HD*
- Nikon Ni Brightfield
- Andor CSU-X1/TIRF*
- Nikon Ts2R Screening

BLSB 914

- Molecular Devices Image Xpress Micro Confocal*

BLSB 1148

- Perkin-Elmer IVIS Lumina

* System has optional **LIVE CELL** steps

NIKON A1R-HD

POWER ON

Power Backup
#1

Press Power button for 1 second until green light begins to flash.
Wait until green light is solid before moving to next step.

Laser Key

Wait 10-15 seconds after.


Confocal
Controller

Power Backup
#2

PC*

Evaluate
Objectives

Software



*If you need just the PC for data retrieval, you must first turn on Power Backup #2

NIKON A1R-HD

LIVE CELL

If you will need incubation for your sample, please turn on the chamber **30 minutes** before you begin.
Incubation controllers are behind the monitor.

Main controller

The settings for each component are set by the factory and will be higher than 37 degrees. Ideal values are printed on the label on the top of the controller. DO NOT ADJUST.

Humidity Bottle

Confirm water level in bottle and refill with Di water if needed.

Remove Stage Adapter

Use the **GREEN** screwdriver stored on the microscope to remove three (3) plastic screws.

Objective Warmer


If using an oil immersion objective.

Install Chamber

Use the **GREEN** screwdriver stored on the microscope to install three (3) plastic screws (hand tight only)

CO2

Turn on tanks, followed by two (2) **RED** valves



If you use incubation – YOU are responsible for returning the system to how you found it – including turning off the CO2.

NIKON A1R-HD

POWER OFF

Clean Objectives

Clean any oil from immersion lenses. If you need additional training on this – contact our staff

Software

Do any necessary processing and close the software

PC


Save any data to removable drives if needed and power down the PC

Power Backup #2

Confocal Controller

Laser Key

Power Backup #1



LEICA STED

POWER ON

Switches 1-3
and key

Wait 15s between switches

Depletion Laser
Key(s)

Only turn on the lasers you
will be using.

Widefield
Power

Evaluate
Objectives

Software



LEICA STED

POWER OFF

Clean Objectives

Clean any oil from immersion lenses. If you need additional training on this – contact our staff

Software

Do any necessary processing and close the software

PC

Save any data to removable drives if needed and power down the PC

Depletion Laser Keys

Key and Switches 1-3



CSU/TIRF

POWER ON

Master Power

Widefield Lamp

PC

Evaluate Objectives

Software



CSU/TIRF

LIVE CELL

If you will need incubation for your sample, please turn on the chamber **30 minutes** before you begin. Incubation controllers are behind the monitor.

Main Power

Install Chamber

CO2

Turn on tanks, followed by two (2) **RED** valves



If you use incubation – YOU are responsible for returning the system to how you found it – including turning off the CO2.

CSU/TIRF

POWER OFF

Clean Objectives

Clean any oil from immersion lenses. If you need additional training on this – contact our staff

Software

PC

Widefield Lamp

Master Power



Widefield Lamp



Master Power



DELTAVISION

POWER ON

Main Power


Bank B Power

PC

Controllers

Objectives

Software



DELTAVISION

LIVE CELL

If you will need incubation for your sample, please turn on the chamber **30 minutes** before you begin. Incubation controllers are behind the monitor.

Control Box

CO2

CO2 tank is housed in BLSB 912.



If you use incubation – YOU are responsible for returning the system to how you found it – including turning off the CO2.

DELTAVISION

Clean Objectives

Clean any oil from immersion lenses. If you need additional training on this – contact our staff

Software

PC

Be sure PC is powered off before turning off other power switches

Controllers

Power Bank B

Main Power

POWER OFF





IMAGE XPRESS



Instrument
Power

PC

Software

Instrument
Power

PC

Software

Sample Heater

CO2 tank

SCREENING

POWER ON

- 1 – System Power
- 2 – Lamp Power
- 3 – Mode Selection

Fluorescence:



- 1 = DAPI
- 2 = GFP
- 3 = RFP
- 4 = Brightfield

USE

- 1 – Choose Filter
- 2 – Adjust Intensity

POWER OFF

- 1 – Lamp Power
- 2 – System Power



Modes:

- DIA = Brightfield/Phase
- EPI = Fluorescence

BRIGHTFIELD

Control Box

Microscope

Camera

Fluorescence*

PC




Control Box

Microscope

Camera

Fluorescence*

PC



IVIS

O2

Pump Switches

Isoflurane

Software



Pump Switches

Software

Isoflurane

Pump Switches

O2

Sanitize

